



Vitamins in Bee Pollen

Ligia Bicudo de Almeida Muradian

Pharmaceutical Science School, University of Sao Paulo, Brazil

Introduction

- Bee pollen → food supplement
- Healthy → high protein content
- Composition → variation according to vegetal species, environment conditions, age and nutritional state of the plant when pollen is under development
- Food Composition Tables / Literature
- Knowledge of chemical composition → specifying the product obtained in different regions
- Method development



Objective

**Study of vitamin composition
of bee pollen loads**



- beta-carotene (provitamin A)
- B Vitamin Complex
- Vitamin C
- Vitamin E

Material and Methods:

Fresh/dried bee pollen (SP/Brazil)
(Storage 0 – 1 year)



- Beta-carotene = Open Column Chromatography
 - B Complex Vitamins
 - Vitamin E
 - Vitamin C (AOAC – titrimetric method)
- }
- HPLC





Methods:

B1

- Conversion from Tiamin to Tiochrome
- Column: C₁₈ (5 µm/125 x 4,0 mm) with pre-column (5 µm/4x4 mm) Lichrospher;
- Detection: fluorescence: Ex 368 nm; Em 440 nm.
- Mobile phase: buffer phosphate pH 7.2 (0,228% KH₂PO₄.3H₂O): dimetilformamide (85:15)
- Flow: 1 mL/min;
- Injection volume: 20 µL;


B2

- Column: C₁₈ (5 µm/125 x 4,0 mm) with pre-column (5 µm/4x4 mm) Lichrospher;
- Detection: fluorescence Ex 450 nm; Em 530 nm.
- Mobile phase: buffer phosphate pH 7.2 (0,228% KH₂PO₄.3H₂O): dimetilformamide (85:15)
- Flow: 1.5 mL/min;
- Injection volume: 20 µL;

Methods:



Vitamin E (alpha-tocopherol)

- Mobile phase: hexane :isopropanol (99:1);
- Flow: 2 mL/min;
- Column: Shim-pack CLC-Sil (M) 25 cm;
- Detection: fluorescence: Ex 295 nm ; Em 330 nm nm.



Chromatogram of alpha-tocopherol

Methods:



Absorption spectrum of beta-carotene

RESULTS

Antioxidant vitamins (dried bee pollen)



Vitamins ($\mu\text{g/g}$)		
C	E	β -carotene
14 \pm 1.62	38.64 \pm 2.72	6.65 \pm 0.29
24 \pm 3.94	16.27 \pm 0.38	3.14 \pm 0.09
40 \pm 11.4	20.54 \pm 1.30	5.05 \pm 0.03
26 \pm 4.88	18.42 \pm 0.47	77.88 \pm 5.01*
27 \pm 2.08	21.06 \pm 1.76	25.19 \pm 1.80
44 \pm 4.49	32.27 \pm 1.53	17.83 \pm 1.91

Results in raw base with mean \pm SD
(n = 3)

* Source of pro-vitamin A
23% RDI for women

From: Melo I. L. P. & Almeida-Muradian, L. B.

Chemical Composition (%)

	Moisture	Proteins	Lipids	Ashes
A	2.56 ± 0.05	24.91 ± 0.88	4.81 ± 0.16	3.17 ± 0.00
B	2.03 ± 0.14	26.45 ± 0.52	4.46 ± 0.20	3.21 ± 0.01
C	2.99 ± 0.04	28.28 ± 0.44	4.53 ± 0.08	3.30 ± 0.03
D	2.18 ± 0.03	21.22 ± 0.18	5.69 ± 0.22	2.90 ± 0.00
E	2.82 ± 0.10	19.98 ± 0.12	5.27 ± 0.32	3.00 ± 0.00
F	1.50 ± 0.05	20.68 ± 0.33	5.18 ± 0.13	2.89 ± 0.01

Results in raw base with mean \pm SD (n = 3)

All results are in accordance with the regulations from; Brazil, France and Argentina

RESULTS

B complex vitamins (dried bee pollen)




	Vitamin B ₁ (mg/25g of pollen)	Vitamin B ₁ (% RDI)	Vitamin B ₂ (mg/25g of pollen)	Vitamin B ₂ (% RDI)
PG	0.25*	20.8	0.62*	47.70
PH	0.17	14.1	0.48*	36.90
PI	0.16	13.3	0.64*	49.20
PJ	0.18*	15.0	0.51*	39.23
PK	0.16	13.3	0.44*	33.85
PL	0.16	13.3	0.44*	33.85
PM	0.20*	16.7	0.51*	39.23
Mean	0.18	15.10	0.52	40.00
RDI	1.2 mg/day		1.3 mg/day	

Results in raw base with mean \pm SD (n = 3)

* Source of the vitamin (B₁ = 17.5 % RDI) (B₂ > 33.85% RDI)

Polinic analysis

(A)
Anadenanthera sp



(B1 and B2)

Eucalyptus sp

(C)

Philodendron sp

(D) type

Areceae

(E) type

Source: MS Dissertation of Karla Oliveira *Myrcia*

Next steps:



- Analyse all vitamins after one year with different storage conditions (dark, light, freezer)
- Analyse other B complex vitamins:
- Vitamers of B₆ (piridoxol, piridoxal and piridoxamine)
- Vitamers of PP (nicotinic acid and nicotinamide)

Acknowledgments

USP Universidade de São Paulo
Brasil



FOOD ANALYSIS LABORATORY – USP - BRAZIL





Thank you!

ligiabi@usp.br